Welch Allyn® Elite Retinoscope Models 18240 and 18245 Welch Allyn® Spot Retinoscope Model 18300

Rx ONLY
REF 773429 DIR 80027580 Ver. A
Revision date: 2021-01

901024 RETINOSCOPE

EC REP and EU IMPORTER
Welch Allyn Limited
Navan Business Park, Dublin Road
Navan, Co. Meath C15 AW22
Ireland

Welch Allyn, Inc.
4341 State Street Road
Skanateles Falls, NY 13153 USA
hillrom.com
Welch Allyn, Inc. is a subsidiary of Hill-Rom Holdings, Inc.

Authorized Australian Sponsor
Welch Allyn Australia Pty. Ltd.
Unit A01, 2-4 Lyndon Park
Macquarie Park, NSW, 2113
Phone 1800 650 083

ENGLISH

Warning and cautions

WARNING The warning statements in this manual identify conditions or practices that could lead to illness, injury, or death. Warning symbols will appear with a grey background in a black and white document.

CAUTION The caution statements in this manual identify conditions or practices that could result in damage to the equipment or other property, or loss of data.

Consult Instructions for use

For information on the origin of these symbols, see the Welch Allyn symbols glossary:
welchallyn.com/symbols/glossary

WARNING Patient injury risk. Use only Welch Allyn approved accessories according to the manufacturer's instructions for use. Using unapproved accessories with the device can affect patient and operator safety and can compromise product performance and accuracy, and void the product warranty.

CAUTION These devices are unlikely to interfere with nearby electronics but may be affected by portable RF communications equipment marked with the symbol. Moving away from the RF transmitter will reduce the interference.

CAUTION These devices may be affected by electrostatic discharge (ESD) events $\geq 15\text{ kV}$. We recommend a use environment with floors made of wood, concrete or other natural materials. If floor coverings are synthetic, relative humidity should be $\geq 30\%$.

CAUTION To minimize lamp-housing temperatures, time should not exceed 2 minutes with off time not less than 10 minutes.

CAUTION Federal Law restricts this device to sale by or order of a physician.

CAUTION While no acute optical radiation hazards have been identified for retinoscopes, it is recommended that the intensity of light directed into the patient level which is necessary for diagnosis. Infants, aphakes, and persons with diseased eyes will be at greater risk. The risk may also be increased if the person being examined has had any exposure with the same instrument or any other ophthalmic instrument using a visible light source during the previous 24 hours. This will apply particularly if the eye has been exposed to retinal photography. The intended use of this device is for routine ophthalmic exams on the order of less than 60 seconds per eye. Significant use of this device beyond its intended use is not recommended; it may cause harm to the eyes.

CAUTION The photochemical weighted retinal irradiances (EB and EA) give a measure of the potential hazard that exists for a beam of light to cause photochemical damage to the retina. EB gives a measure for eyes in which a crystalline lens is in place. EA applies to eyes in which the crystalline lens has been removed (aphakes) and has not been replaced by a UV blocking lens or for the eyes of very young children. The photochemical retinal exposure dose limit is 2.54 Joules/cm² based on calculations arising from the American Conference of Governmental Industrial Hygienists (ACGIH) Threshold Limit Values (TLV) limit divided by the photochemically weighted spectral irradiances gives the maximum permissible exposure (m.p.e.).

m.p.e. = 434 seconds at 40 cm
268 seconds at 66 cm
m.p.e. (Aphake) = 267 seconds at 40 cm
164 seconds at 66 cm

Transport/storage/operating

-20°C - 49°C, 95% RH Max

Repair and maintenance

No disinfection is required before service.


Standards and compliance

The device complies with the following standards:

ISO 12865, ISO 15004-1, ISO 10993-1, IEC 60601-1 and IEC 60601-1-2.

Country-specific standards are included in the applicable Declaration of Conformity.

The spectrum as it exits the device is shown below.

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