Hillrom

WELCH ALLYN Elite Retinoscope
18240, 18245, 18300

Instructions for use

REF 775109, 80029246 Ver. B, Revision date: 2022-11

901024 RETINOSCOPE Rx ONLY

Baxter, Hillrom, and Welch Allyn are trademarks of Baxter International, Inc., or its subsidiaries.

Welch Allyn, Inc. ("Welch Allyn") assumes no responsibility for any injury to anyone that may result from (i) failure to properly use the product in accordance with the instructions, cautions, warnings, or statement of intended use published in this manual, or (ii) any illegal or improper use of the product.

Welch Allyn, Inc.
4341 State Street Road
Schenectady, NY 12151 USA
hillrom.com

Welch Allyn, Inc. is a subsidiary of Hill-Rom Holdings, Inc.

EC REP and EU IMPORTER

Welch Allyn Limited
Navan, Co. Meath C15 AW22
Ireland

Welch Allyn Australia Sponsor

1 Buxton Drive,
Old Toongabbie NSW 2146
Australia

ENGLISH

Instructions for use

Thank you for purchasing the WELCH ALLYN 3.5V Retinoscope. This instrument has been designed to meet the needs of today's practitioners and incorporates features not found on any other retinoscopy.

1. External Fixing Sleeve - Unique planetary gear system allows for easy adjustment no matter what size hand or how instrument is held. Continuous 360° rotations. Maintains the same plane of focusing during rotation.

2. Improved Spot Optics - Outer halogen lamp provides 50% more intense yellow-green lamps. The reflex is now easier and easier to see in all patient's eyes. Retinoscopy can be done faster and more accurately.

3. Dust-free Optics - New housings and glass cover on the front keep the instrument cleaner longer.

4. Crossed linear Polarizing Filter - Dramatically reduces glare from lenses. Allows retinoscopy to be performed closer to the axis of the correcting lenses.

5. Fixation Cards - New cards that easily attach increase the ease with which dynamic retinoscopy is performed.

6. Improved Optics - Glare and shadows have been eliminated for a clearer and more precise view.

7. Interchangeability - By simply changing the lamp, the streak retinoscopy can be converted to a spot retinoscopy.

Intended use

A retinoscope is an AC-powered or battery-powered device intended to measure the refraction of the eye by illuminating the retina and noting the direction of movement of the light on the retina surface and note the refraction by the eye of the emergent eye. Works in conjunction with your WELCH ALLYN Retinoscopes are intended for use by clinicians and medically qualified personnel. WELCH ALLYN Retoscopes are intended for use in both ambulatory and acute care settings.

Indications for use

The Retinoscope is intended to aid in measuring the refraction of the eye.

Symbol descriptions

For information on the origin of these symbols, see the Welch Allyn symbols glossary: welchallyn.com/symbols/glossary.

MD	Medical device	#	Product Identifier
	Manufacturer		Authorized Representative in the European Community
	Lot code		Consult Instructions for use
	Reorder number		Global Trade Item Number
	Temperature limitation		Humidity limitation

Rx ONLY	Prescription only or "For use by or on the order of a licensed medical professional"		Recyclable
	Separate collection of Electrical and Electronic Equipment. Do not dispose as unsorted municipal waste.		

CAUTION The caution statements in this manual identify conditions or practices that could result in damage to the equipment or other property, or loss of data.

WARNING The warning statements in this manual identify conditions or practices that could lead to illness, injury, or death. Warning symbols will appear with a gray background in a black and white document.

Residual risk

This product complies with relevant electro-magnetic interference, mechanical safety, performance, and biocompatibility standards. However, the product cannot completely eliminate potential patient or user harm from the following:

- Harm or device damage associated with electro-magnetic hazards,
- Harm from mechanical hazards,
- Harm from device, function, or parameter unavailability,
- Harm from misuse error, such as inadequate cleaning, and/or
- Harm from device exposure to biological triggers that may result in a severe systemic allergic reaction.

Product use

The operator will note that the width of the streak varies as the sleeve is raised and lowered (see Figure 1). When the operating sleeve is in the lowest position the light rays emitted are slightly divergent. However, divergent rays that will recombine to focus. As the sleeve is raised, the streak focuses. With the sleeve all the way up, the retinoscope acts with a concave mirror effect, where the light rays cross and then diverge. Because the rays cross, the eye's reflex moves in opposite directions with the concave mirror effect as compared to the plano mirror effect.

Throughout this IFU, we will use the plano mirror effect unless specified.

The rotary movement of the control sleeve mechanism allows the streak to rotate 360° to ascertain the axis of astigmatism (Figure 1).

Preliminary steps

- Set the sleeve in its lowest position (plano mirror effect).
- Position the sleeve 2/3 meters (26') from the patient. This distance is a working distance of 1.50 (compared as the reciprocal of working distance in meters). Working distance and motion may be varied to suit the practitioner's needs. (In this instruction book, the 2/3-meter (26') working distance is assumed. Different working distances can be used, but remember to adjust for your patient's distance.)
- With the refracting instrument in place, direct the patient's attention to a fixation spot at 1.5 feet or more from the eye and align the streak vertically.

4. Observe the "reflex" which will appear as in Figure 2, providing no oblique astigmatism is present. If oblique astigmatism is present, the reflex will appear more like Figure 3, where the reflex does not appear vertical.

5. Move the vertical streak horizontally across the pupil and back again and observe whether the reflex moves in the same direction as the streak or in the opposite direction.

Figure 2 Figure 3

6. Rotate the control sleeve until the streak is horizontal and move the streak vertically. The reflex will appear as in Figures 4 or 5.

Figure 4 Figure 5

7. If the streak and the reflex move in the same direction with no lens in the refractive apparatus, the refraction is one of these:

- Hyperopia;
- Emmetropia;
- Myopia of less than 1.50 diopters.

If the reflex moves in the opposite direction, the error is myopia greater than 1.50 diopters.

Deferring refractive error by neutralization

Before starting, make sure the eye not being refracted has some "against" motion using the plano mirror effect. This will blur vision to prevent accommodation. This is true for neutral refraction is noticed initially, place about a +1.00 sphere before the eye once neutral motion is seen.

Neutralizing with spheres only

1. Change sphere in the minus direction until the reflexes in all axes have "with" motion.

2. Adjust in the plus direction until the reflex fills the pupil in one meridian and all motion has stopped. This will be one of the principal meridians if astigmatism is present. That meridian is then said to be neutralized.

3. Test for neutralization by one of these methods:

- Move the sleeve all the way up (converge mirror position); the reflex should also appear neutralized;
- Move closer to the patient and "with" motion should return; move away and "against" motion should appear;
- Place an extra +0.25 sphere in the apparatus and "against" motion should appear;

4. Repeat the neutralization in the meridian 90° away.

Figure 6

Locating the axis of astigmatism
Two phenomena occur when determining the axis of astigmatism: break and streak. A streak is observed when the streak is aligned with a principal meridian of the astigmatism (Figures 3 and 5). The streak will be aligned with a principal meridian when the break effect disappears and the width of the reflex is narrowest (and it appears brightest, Figure 6).

This LIMITED WARRANTY IS PROVIDED INSTEAD OF ANY OTHER IMPLIED PRODUCT QUALITY WARRANTIES AND IT IS THE SOLE WARRANTY OF THE PRODUCTS.

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