



HILL-ROM® 900 BED

Improve safety. Simplify tasks.

FEATURES AND BENEFITS

The lowering of the split siderails can be activated using one hand and feature a smooth, rapid self-lowering mechanism allowing for quick access to patients

Stable support is provided during bed exit via an ergonomically-placed side egress grip handle

Intuitive bed controls require minimal in-service staff training and help minimise the risk of user's error

One-Touch Side Egress positions the bed for safer and easier bed exit

AutoContour™ and SlideGuard™ minimise patient migration in bed, reducing abdominal compression during bed articulation

EasyChair™ and fl at position function offers patients the ability to quickly reach a comfortable seated position

Low height indicator warns the caregiver if the bed is not in the lowest position, contributing to a safer patient environment

5th castor facilitates easier bed transport in any direction

One-Touch Shock positioning fl attens the upper deck and places the bed into a Trendelenburg position, should patient complications occur

Intelligent night light enhances operational safety for both patients and caregivers

Alerts caregivers of potential bed exit via the 3-Mode Bed Exit Alarm, available with remote monitoring capabilities*

* Option available in selected markets.

TECHNICAL SPECIFICATIONS

Safe Working Load 220 kg 250 kg	+ o
Bed Weight	144 kg
Height Adjustment	38.9 to 76.8 cm
SlideGuard Sliding Backrest	up to 12 cm
Protection Class	IPX4
Bed Angles Head Section Thigh Section Foot Section Trendelenburg/Reverse Trendelenburg	0 to 65° 0 to 28° -22 to -3° -17 to +17°
Sleep Deck Width Length	90 cm 200 cm
Sleep Deck Surface	High Pressure Laminate (HPL)
Siderail Height	39.3 cm
Castors 150 mm Integral 125 mm Double Band** 150 mm Double Band	+ o o
Steering Castor At Head End At Foot End	+ o
5th Castor	o
Brake and Steer Pedals At Foot End At Foot and Head Ends**	+ o
Brake-Off Indicator	o
3-Mode Bed Exit Alarm	o
Nurse Call	o
Controls Integrated in Siderails Not Integrated in Siderails (Pendant) Pendant** Control on Flexible Arm** Caregiver Foot End Pendant** Bilateral High-Low Foot Control	+ o o o o o
Head Section Mobile Stationary	+ o
Low Height Indicator**	+
EasyChair and Flat Position Function	+
Intelligent Night Light	o
One-Touch Shock Positioning**	o
Bed Extension with Linen Holder	o

+ = Included; O = Optional

** Only available in combination with selected options. For more information please contact your local Hillrom representative.



Optional intelligent night light



Low height indicator and One-Touch Side Egress positioning



Optional 125 mm double band castors**



3-Mode Bed Exit Alarm



Care grip handles and integrated controls



Optional bilateral high-low foot control with automatic lock-out and 5th steering castor

ADDITIONAL STANDARD FEATURES

- Bilateral Mechanical CPR Lever
- 6 Attachment Points for Patient Restraints
- Lifting/IV Pole Holders
- Bilateral Accessory Holder
- Electric Cord Holder
- Adjustable Mattress Retainers
- Battery
- Removable Endboards
- Line-of-Site Backrest and Trendelenburg Indicators
- One-Touch Side Egress Positioning
- Central Brake
- Service Required Indicator

ABOUT HILLROM

Hillrom is a global medical technology leader whose 10,000 employees have a single purpose: enhancing outcomes for patients and their caregivers by advancing connected care. Around the world, our innovations touch over 7 million patients each day. They help enable earlier diagnosis and treatment, optimize surgical efficiency and accelerate patient recovery while simplifying clinical communication and shifting care closer to home. We make these outcomes possible through connected smart beds, patient lifts, patient assessment and monitoring technologies, caregiver collaboration tools, respiratory care devices, advanced operating room equipment and more, delivering actionable, real-time insights at the point of care. Learn more at hillrom.com.



Hillrom[™]

hillrom.com

The Hill Rom 900 bed is intended for use in intensive, acute and ambulatory care settings for patients with a weight ≥ 40 kg, height ≥ 146 cm and BMI ≥ 17 .
Class I
Manufacturer: Hill-Rom S.A.S. – BP 14 ZI du Talhouët 56330 Pluvigner France

This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.