



The Welch Allyn Green Series™ 777 Wall Transformer

Helps reduce costs with an energy-efficient product.

- Energy Star certified 5 W power supply
- Helps reduce the risk of cross-contamination with no on/off switch
- Compatible with all 3.5 V Welch Allyn instrument heads
- Backwards compatible to Welch Allyn 767 Wall Transformers and Integrated Wall Systems

Specifications

| | | |
|---|--|---|
| Dimensions | Width: 180 mm | |
| | Height: 148 mm | |
| Environment | Depth: 111 mm | |
| | Weight: 0.76 kg | |
| Safety Class | Power Handle Cord Length (extended): 3.66 m | |
| | Operating: | Temperature limitation: 15°C to 40°C Humidity limitation: 15% to 95% Atmospheric limitation: 700 hPa to 1060 hPa |
| Electrical Power supply | Transport/storage: | Temperature limitation: -20°C to 50°C Humidity limitation: 15% to 95% Atmospheric limitation: 700 hPa to 1060 hPa |
| | II Intermittent Operation | |
| Materials | Input: 100–240 VAC, 50-60 Hz, 0.18 A | Output: 5 VDC 1 A |
| Approvals | ABS and elastomer rubber. Latex free. | |
| Warranty | IEC 60601-1, UL 60601-1, IEC 60601-1-2, CAN/CSA C22.2 No 601.1-M90, AS/NZS3200 | |
| Guidance and manufacturer's declaration | 10-year warranty | |
| | The guidance and manufacturer's declaration information is available at www.welchallyn.com/GreenSeries777WallTransformer | |

Accessories

| | |
|-----------|------------------------------|
| 713262 | 5 W PLATFORM POWER SUPPLY |
| 713264 | 5 W ADAPTER PLUG UK |
| #104714 | 0.457 m Twin USB Power Cable |
| #76700 | Locking Collar |
| #76710-12 | Mobile Stand |

Welch Allyn UK Ltd. Cublington Road, Buckinghamshire, HP22 4ND
Telephone: 0207 365 6780 | Fax: 0207 365 9694

United Kingdom: +44 207 365 6780

For complete details on the Welch Allyn Green Series Wall Transformer and Integrated Wall System, visit us online at www.welchallyn.com

