



Welch Allyn KleenSpec® Disposable Vaginal Specula and Cordless Illumination System



An advancement in women's health technology
that's nothing short of brilliant



WelchAllyn®

Welch Allyn KleenSpec® Disposable Vaginal Specula
and Cordless Illumination System—
now in every size



590 Series KleenSpec® Disposable Vaginal Specula

Enhance Patient Comfort

- The extra small speculum allows you to perform a comfortable exam with virginal/pediatric patients as well as post-menopausal and post-hysterectomy women—no more referrals or uncomfortable exams using metal specula or the wrong size
- Smooth, molded acrylic designed to prevent pinching or binding of tissue for greater patient comfort

Help Reduce Cross-Contamination

- Single-use speculum helps to reduce cross-contamination for improved patient safety while eliminating the time and expense of reprocessing metal specs

Improve Staff Productivity

- Features a wider handle design for improved maneuverability and better balance during exams
- A strong grip and a yoke that's more flexible than low-cost substitutes



Look for the distinctive Welch Allyn blue, green, white and yellow colored handles

“ We have been using the Welch Allyn cordless Illumination system with disposable vaginal speculums for all of our patient exams since 2008 and couldn't be happier. This light is brighter than any I have ever seen and shows true tissue color. If you forget to charge the light, it is easy to grab a replacement from another room or use a spare. The light is small, easy and quick to use which truly makes us more efficient. I would recommend this system to any clinicians looking to improve the quality of their care.

Shelley Kolton, MD / Downtown Women OB/Gyn Associates / New York City, NY



KleenSpec Cordless Illumination System

- LED technology as opposed to halogen, for a whiter, brighter light and true tissue color during an exam
- No bulb replacement necessary for a reduced cost of ownership
- Features a cordless design so you'll never have to worry about cords getting in your way, breaking or needing to be cleaned, speeding the exam and giving you a full range of motion
- Push-button on/off; the illumination system easily fits into every size 590 Series speculum
- Newly designed light pipe provides a homogeneous spot and eliminates glare back so all of the light is projected forward
- Uniform spot—no dark or hot spots—for enhanced visualization of the exam area
- Produces much less heat than traditional lamps and has a continuous on time of 80 minutes
- A lithium ion battery provides a recharge time of five hours



KleenSpec Vaginal Specula—Premium 590 Series

- 590XS** KleenSpec 590 Series Premium Disposable Vaginal Specula, **X-Small (24/box; 4 boxes/case)**
- 59000** KleenSpec 590 Series Premium Disposable Vaginal Specula, **Small (24/box; 4 boxes/case)**
- 59001** KleenSpec 590 Series Premium Disposable Vaginal Specula, **Medium (24/box; 4 boxes/case)**
- 59004** KleenSpec 590 Series Premium Disposable Vaginal Specula, **Large (18/box; 4 boxes/case)**

KleenSpec Vaginal Specula—Premium 590 Series with Smoke Tube

- 59005** KleenSpec 590 Series Premium Disposable Vaginal Specula with Smoke Tube, **Small (12/box; 4 boxes/case)**
- 59006** KleenSpec 590 Series Premium Disposable Vaginal Specula with Smoke Tube, **Medium (12/box; 4 boxes/case)**

KleenSpec Cordless Illumination System

- 80015** Complete KleenSpec Cordless Illuminator with Charging Station, EU and UK plugs (conforms to IEC 60601-1)
- 80000** KleenSpec Cordless Illuminator
- 74015** Charging Station for Cordless Illuminator, EU and UK plugs
- 59010** Disposable Sheath for Cordless Illumination System (**100/box; 5 boxes/case**)

KleenSpec Vaginal Specula—Premium 580 Series w/Sheath

- 58000S** KleenSpec 580 Series Premium Disposable Vaginal Specula with Sheath, **Small (18/box; 4 boxes/case)**
- 58001S** KleenSpec 580 Series Premium Disposable Vaginal Specula with Sheath, **Medium (18/box; 4 boxes/case)**
- 58004S** KleenSpec 580 Series Premium Disposable Vaginal Specula with Sheath, **Large (12/box; 4 boxes/case)**

Contact your Welch Allyn representative today or visit <http://emeai.welchallyn.com> to learn more

Welch Allyn EMEA Customer Care
IDA Business Park - Johnstown
Navan, County Meath, Ireland
(p) +353 46 906 7700
europe@welchallyn.com
EMEAI.WELCHALLYN.COM



WelchAllyn[®]

Hill-Rom[®]