

October 9, 2023

Dear Valued Customer,

We are pleased to announce that Baxter is releasing a new barcode scanner compatible with the **Welch Allyn Connex** Spot Monitor. The new barcode scanner is the “**Zebra CS6080** Scanner and Holster Kit”, part number 7000-CS6080. The 7000-CS6080 barcode scanner can read 1D and 2D barcodes; it does not have RFID functionality.

The 7000-CS6080 scanner is only compatible with the **Connex** Spot Monitor at this time and requires **Connex** Spot Monitor Firmware 1.40.10 A0002 or greater for use. This will be an additional barcode scanner offering, Baxter is not discontinuing the sale of any of its other current scanner offerings.

The 7000-CS6080 barcode scanner will ship with a coiled cord and barcode scanner holster, but replacement coiled cords and barcode scanner holsters will also be available for individual purchase.

New Barcode Scanner

Zebra CS6080 Scanner and Holster Kit, Part Number: 7000-CS6080



New Barcode Scanner Replacement Parts

Universal Scanner Holster Kit, Part Number: 108781

- Note: Part 108781 is compatible with the 7000-CS6080, 7000-916HS, 7000-916HSR, 6000-916HS, and 6000-916HSR Scanners

Zebra CS6080 Coiled Cord Replacement, Part Number: 108782

- Note: Part 108782 is currently only compatible with the 7000-CS6080 Scanner

If you desire the 7000-CS6080 barcode scanner to be configured, please contact our customer service department.

To obtain more information about the 7000-CS6080 barcode scanner, please visit Hillrom.com or contact your local Baxter sales representative or customer care for additional details.

Sincerely,



James Bennett
Associate Marketing Manager, Global Marketing
Patient Monitoring
Baxter Healthcare Corporation

Baxter, Connex and Welch Allyn are trademarks of Baxter International Inc. or its subsidiaries.
Zebra is a trademark of Zebra Technologies Corp.